


Smooth Round Low Profile Gel Implant

Cohesive I™



Cohesive I™ gel is the standard gel used in Mentor's implants. This is Mentor's softest gel and has been preferred for years in Round Low, Moderate, Moderate Plus, High Profile and Becker implants.

- Smooth Surface, available in Low Profile
- Offered in a Wide Range of Sizes
- Features Mentor's Proprietary Cohesive Gel Formulation
- Shell Barrier Layer Designed to Reduce Gel Bleed
- Lifetime Product Replacement Guarantee

CE Not for distribution in the U.S.



the power to transform®

MENTOR




Smooth Round Low Profile Gel Breast Implant

Smooth Round Low Profile, Cohesive I™

Vol. cc	Diam. cm	Proj. cm	Catalog Number	Resterilizable Gel Sizer
100	9,4	1,9	310JSP1100	RSZ-1100
125	10,1	2,0	310JSP1125	RSZ-1125
150	10,7	2,2	310JSP1150	RSZ-1150
175	11,3	2,3	310JSP1175	RSZ-1175
200	11,8	2,4	310JSP1200	RSZ-1200
225	12,3	2,5	310JSP1225	RSZ-1225
250	12,7	2,6	310JSP1250	RSZ-1250
275	13,1	2,7	310JSP1275	RSZ-1275
300	13,5	2,7	310JSP1300	RSZ-1300
325	13,9	2,8	310JSP1325	RSZ-1325
350	14,2	2,9	310JSP1350	RSZ-1350
375	14,5	2,9	310JSP1375	RSZ-1375
400	14,8	3,0	310JSP1400	RSZ-1400
450	15,4	3,1	310JSP1450	RSZ-1450
500	16,0	3,3	310JSP1500	RSZ-1500
550	16,5	3,4	310JSP1550	RSZ-1550
600	17,0	3,5	310JSP1600	RSZ-1600

Note: Individual implant dimensions may vary slightly in products of this type. Not all units will conform exactly to the dimensions noted above.



INTERNATIONAL OPERATIONS



201 Mentor Drive
Santa Barbara, CA 93111 USA
+1 (805) 879-6000
Fax: +1 (805) 967-7108
www.mentorcorp.com

Manufacturer:
Mentor Medical Systems B.V.
Zernikedreef 2
2333 CL Leiden
The Netherlands

Mentor Medical Systems, Australia, Pty. Ltd. Tel: +61 2 9251 8111 Toll Free: +1 800 338 160 Fax: +61 2 9251 8333
Mentor Medical Systems, Canada, Inc. Tel: +1 905 725 7763 Toll Free: +1 800 668 6069 Fax: +1 905 725 7340
Mentor Benelux B.V. Belgium: Tel: +32 15 45 91 56 Fax: +32 15 45 91 57
Mentor Benelux B.V. The Netherlands: Tel: +31 71 524 99 00 Fax: +31 71 524 99 48
Mentor Deutschland GmbH Tel: +49 811 600 50 0 Fax: +49 811 600 50 11
Mentor Medical Systems France, S.A. Tel: +33 1 46 01 30 54 Fax: +33 1 46 01 30 55
Mentor Medical Systems Iberica, S.L. Tel: +34 91 562 2700 Fax: +34 91 728 3367
Mentor Medical Italia, S.r.l. Tel: +39 0 2 880 7761 Fax: +39 0 2 7209 3710
Mentor Medical Systems Ltd., U.K. Tel: +44 1635 511800 Fax: +44 1635 511810