


Aseptic Transfer System



The Aseptic Transfer System is indicated for the aseptic transfer of sterile, isotonic saline for injection from an I.V. container into an inflatable breast implant or tissue expander.

Aseptic Transfer System



Carton contains four sterile kits. Each kit contains the following:

Quantity	Catalog Number	Description
1	350-8400	Aseptic Transfer System
		60 cc syringe (sterile, disposable)
		Piercing spike
		48" tubing (saline container to syringe)
		Dual check valve
		10" tubing (syringe to luer connector)
		Luer connector (attaches to implant fill-tube)



201 Mentor Drive
Santa Barbara, CA 93111 USA

(805) 879-6000

Customer Service:
(800) 235-5731
Fax: (805) 967-7108

www.mentorcorp.com

Manufacturer:
Mentor
3041 Skyway Circle North
Irving, TX 75038 USA

International Operations

Mentor Medical Systems, Australia, Pty. Ltd. Tel: +61 2 9251 8111 Toll Free: +1 800 338 160 Fax: +61 2 9251 8333

Mentor Medical Systems, Canada, Inc. Tel: +1 905 725 7763 Toll Free: +1 800 668 6069 Fax: +1 905 725 7340

Mentor Deutschland GmbH Tel: +49 811 600 50 0 Fax: +49 811 600 50 11

Mentor Medical Systems France, S.A. Tel: +33 1 30 07 07 60 Fax: +33 1 30 07 07 67

Mentor Medical Systems Iberica, S.L. Tel: +34 91 562 2700 Fax: +34 91 728 3367

Mentor Medical Italia, S.r.l. Tel: +39 0 2 880 7761 Fax: +39 0 2 7209 3710

Mentor Medical Systems Ltd., U.K. Tel: +44 1635 511800 Fax: +44 1635 511810