


# Contour Profile® Gel Family – The 300 Series

Cohesive III™



Filled with Mentor's most cohesive gel, providing shape retention with a pleasing level of firmness for optimal aesthetic results.

- Three Different Heights to Match Every Body Type
- Three Different Projections to Maximize Post-Operative Possibilities
- Features Mentor's Proprietary Cohesive Gel Formulation
- Aggressive Siltex™ Texturing
- Orientation Marks for Easier Surgical Placement

Not for distribution in the U.S.



# Contour Profile® Gel Family – The 300 Series



CPG™ NAMING LEGEND		
Cohesive Level	Implant Height	Implant Projection
III 3	Tall 3	High 3
	Medium 2	Moderate <i>Plus</i> 2
	Low 1	Moderate 1
KEY: Digit 1 = Cohesive Level Digit 2 = Height of Implant Digit 3 = Projection of Implant		

## CPG™ 332 Cohesive III™, Tall Height, Moderate *Plus* Profile

Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
145 cc	9,0 cm	9,4 cm	3,8 cm	334-0909	RSZ-0909
175 cc	9,5 cm	9,9 cm	4,0 cm	334-0959	RSZ-0959
205 cc	10,0 cm	10,4 cm	4,2 cm	334-1009	RSZ-1009
235 cc	10,5 cm	10,9 cm	4,4 cm	334-1059	RSZ-1059
270 cc	11,0 cm	11,5 cm	4,7 cm	334-1109	RSZ-1109
305 cc	11,5 cm	12,0 cm	4,9 cm	334-1159	RSZ-1159
350 cc	12,0 cm	12,5 cm	5,1 cm	334-1209	RSZ-1209
395 cc	12,5 cm	13,0 cm	5,3 cm	334-1259	RSZ-1259
445 cc	13,0 cm	13,5 cm	5,5 cm	334-1309	RSZ-1309
495 cc	13,5 cm	14,1 cm	5,7 cm	334-1359	RSZ-1359
555 cc	14,0 cm	14,6 cm	5,9 cm	334-1409	RSZ-1409
615 cc	14,5 cm	15,1 cm	6,1 cm	334-1459	RSZ-1459
680 cc	15,0 cm	15,6 cm	6,3 cm	334-1509	RSZ-1509

## CPG™ 333 Cohesive III™, Tall Height, High Profile

Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
180 cc	9,0 cm	9,4 cm	4,6 cm	334-0904	RSZ-4090
215 cc	9,5 cm	9,9 cm	4,8 cm	334-0954	RSZ-4095
250 cc	10,0 cm	10,4 cm	5,1 cm	334-1004	RSZ-4100
290 cc	10,5 cm	10,9 cm	5,3 cm	334-1054	RSZ-4105
330 cc	11,0 cm	11,5 cm	5,6 cm	334-1104	RSZ-4110
380 cc	11,5 cm	12,0 cm	5,8 cm	334-1154	RSZ-4115
430 cc	12,0 cm	12,5 cm	6,0 cm	334-1204	RSZ-4120
485 cc	12,5 cm	13,0 cm	6,2 cm	334-1254	RSZ-4125
545 cc	13,0 cm	13,5 cm	6,5 cm	334-1304	RSZ-4130
610 cc	13,5 cm	14,1 cm	6,7 cm	334-1354	RSZ-4135
680 cc	14,0 cm	14,6 cm	6,9 cm	334-1404	RSZ-4140
755 cc	14,5 cm	15,1 cm	7,1 cm	334-1454	RSZ-4145

## CPG™ 321 Cohesive III™, Medium Height, Moderate Profile

Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
120 cc	9,0 cm	8,5 cm	3,3 cm	354-0908	RSZ-0908
135 cc	9,5 cm	8,9 cm	3,5 cm	354-0958	RSZ-0958
155 cc	10,0 cm	9,4 cm	3,7 cm	354-1008	RSZ-1008
180 cc	10,5 cm	9,9 cm	3,8 cm	354-1058	RSZ-1058
215 cc	11,0 cm	10,3 cm	3,9 cm	354-1108	RSZ-1108
245 cc	11,5 cm	10,8 cm	4,0 cm	354-1158	RSZ-1158
280 cc	12,0 cm	11,3 cm	4,2 cm	354-1208	RSZ-1208
315 cc	12,5 cm	11,8 cm	4,4 cm	354-1258	RSZ-1258
355 cc	13,0 cm	12,2 cm	4,6 cm	354-1308	RSZ-1308
395 cc	13,5 cm	12,7 cm	4,7 cm	354-1358	RSZ-1358
440 cc	14,0 cm	13,2 cm	4,9 cm	354-1408	RSZ-1408
480 cc	14,5 cm	13,6 cm	5,0 cm	354-1458	RSZ-1458
530 cc	15,0 cm	14,1 cm	5,2 cm	354-1508	RSZ-1508
640 cc	16,0 cm	15,0 cm	5,6 cm	354-1608	RSZ-1608
775 cc	17,0 cm	16,0 cm	5,9 cm	354-1708	RSZ-1708

## CPG™ 322 Cohesive III™, Medium Height, Moderate *Plus* Profile

Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
140 cc	9,0 cm	8,5 cm	3,8 cm	334-0905	RSZ-0905
165 cc	9,5 cm	8,9 cm	4,0 cm	334-0955	RSZ-0955
195 cc	10,0 cm	9,4 cm	4,2 cm	334-1005	RSZ-1005
225 cc	10,5 cm	9,9 cm	4,4 cm	334-1055	RSZ-1055
255 cc	11,0 cm	10,3 cm	4,7 cm	334-1105	RSZ-1105
295 cc	11,5 cm	10,8 cm	4,9 cm	334-1155	RSZ-1155
330 cc	12,0 cm	11,3 cm	5,1 cm	334-1205	RSZ-1205
375 cc	12,5 cm	11,8 cm	5,3 cm	334-1255	RSZ-1255
420 cc	13,0 cm	12,2 cm	5,5 cm	334-1305	RSZ-1305
475 cc	13,5 cm	12,7 cm	5,7 cm	334-1355	RSZ-1355
525 cc	14,0 cm	13,2 cm	5,9 cm	334-1405	RSZ-1405
585 cc	14,5 cm	13,6 cm	6,1 cm	334-1455	RSZ-1455
650 cc	15,0 cm	14,1 cm	6,3 cm	334-1505	RSZ-1505

## CPG™ 323 Cohesive III™, Medium Height, High Profile

Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
165 cc	9,0 cm	8,5 cm	4,6 cm	334-0902	RSZ-0902
195 cc	9,5 cm	8,9 cm	4,8 cm	334-0952	RSZ-0952
225 cc	10,0 cm	9,4 cm	5,1 cm	334-1002	RSZ-1002
260 cc	10,5 cm	9,9 cm	5,3 cm	334-1052	RSZ-1052
300 cc	11,0 cm	10,3 cm	5,6 cm	334-1102	RSZ-1102
345 cc	11,5 cm	10,8 cm	5,8 cm	334-1152	RSZ-1152
390 cc	12,0 cm	11,3 cm	6,0 cm	334-1202	RSZ-1202
440 cc	12,5 cm	11,8 cm	6,2 cm	334-1252	RSZ-1252
495 cc	13,0 cm	12,2 cm	6,5 cm	334-1302	RSZ-1302
555 cc	13,5 cm	12,7 cm	6,7 cm	334-1352	RSZ-1352
620 cc	14,0 cm	13,2 cm	6,9 cm	334-1402	RSZ-1402
685 cc	14,5 cm	13,6 cm	7,1 cm	334-1452	RSZ-1452

## CPG™ 312 Cohesive III™, Low Height, Moderate *Plus* Profile

Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
125 cc	9,0 cm	8,0 cm	3,8 cm	334-0907	RSZ-0907
145 cc	9,5 cm	8,4 cm	4,0 cm	334-0957	RSZ-0957
170 cc	10,0 cm	8,8 cm	4,2 cm	334-1007	RSZ-1007
195 cc	10,5 cm	9,3 cm	4,4 cm	334-1057	RSZ-1057
225 cc	11,0 cm	9,7 cm	4,7 cm	334-1107	RSZ-1107
255 cc	11,5 cm	10,2 cm	4,9 cm	334-1157	RSZ-1157
290 cc	12,0 cm	10,6 cm	5,1 cm	334-1207	RSZ-1207
330 cc	12,5 cm	11,1 cm	5,3 cm	334-1257	RSZ-1257
370 cc	13,0 cm	11,5 cm	5,5 cm	334-1307	RSZ-1307
415 cc	13,5 cm	12,0 cm	5,7 cm	334-1357	RSZ-1357
465 cc	14,0 cm	12,4 cm	5,9 cm	334-1407	RSZ-1407
515 cc	14,5 cm	12,8 cm	6,1 cm	334-1457	RSZ-1457
570 cc	15,0 cm	13,3 cm	6,3 cm	334-1507	RSZ-1507
690 cc	16,0 cm	14,2 cm	6,8 cm	334-1607	RSZ-1607

Note: Individual implant dimensions may vary slightly in products of this type. Not all units will conform exactly to the dimensions noted above.

\*Resterilizable up to 10 times.

### INTERNATIONAL OPERATIONS

**Mentor Medical Systems, Australia, Pty. Ltd.** Tel: +61 2 9251 8111 Toll Free: +1 800 338 160 Fax: +61 2 9251 8333  
**Mentor Medical Systems, Canada, Inc.** Tel: +1 905 725 7763 Toll Free: +1 800 668 6069 Fax: +1 905 725 7340  
**Mentor Deutschland GmbH** Tel: +49 811 600 50 0 Fax: +49 811 600 50 11  
**Mentor Medical Systems Iberica, S.L.** Tel: +34 91 562 2700 Fax: +34 91 728 3367  
**Mentor Medical Italia, S.r.l.** Tel: +39 0 2 880 7761 Fax: +39 0 2 7209 3710  
**Mentor Medical Systems Ltd., U.K.** Tel: +44 1635 511800 Fax: +44 1635 511810  
**Perouse Plastie, S.A.S.** Tel: +33 3 44 08 26 26 Fax: +33 3 44 08 26 01



201 Mentor Drive  
 Santa Barbara, CA 93111 USA  
 Tel: +1 805 879 6000  
 Fax: +1 805 967 7108  
 www.mentorcorp.com

Manufacturer:  
 Mentor Medical Systems B.V.  
 Zernikedreef 2  
 2333 CL Leiden  
 The Netherlands