

Sterile Resterilizable* Gel Breast Implant Sizers



• Resterilizable Gel Breast Implant Sizers for Volume Determination

• Shipped Sterile – *Resterilizable up to 10 Times

• Enhanced Autoclaving Protocols

• Sizer Clearly Marked for Simple Indication of the Volume Size and the Style

• For Temporary Intraoperative Use

• Dimensions Compatible with our [MemoryGel®](#) Product Lines



Not for distribution in the U.S.



Sterile Resterilizable* Gel Breast Implant Sizers

CPG™ 331^{NEW} Cohesive III™, Tall Height, Moderate Projection

Vol.	Width	Height	Proj.	Catalog #
125 cc	9,0 cm	9,2 cm	3,2 cm	RSZ-0903
150 cc	9,5 cm	9,7 cm	3,3 cm	RSZ-0953
175 cc	10,0 cm	10,2 cm	3,4 cm	RSZ-1003
200 cc	10,5 cm	10,7 cm	3,5 cm	RSZ-1053
230 cc	11,0 cm	11,3 cm	3,6 cm	RSZ-1103
265 cc	11,5 cm	11,8 cm	3,7 cm	RSZ-1153
300 cc	12,0 cm	12,3 cm	3,9 cm	RSZ-1203
340 cc	12,5 cm	12,8 cm	4,0 cm	RSZ-1253
380 cc	13,0 cm	13,3 cm	4,1 cm	RSZ-1303
425 cc	13,5 cm	13,8 cm	4,3 cm	RSZ-1353
475 cc	14,0 cm	14,3 cm	4,5 cm	RSZ-1403
530 cc	14,5 cm	14,8 cm	4,6 cm	RSZ-1453
585 cc	15,0 cm	15,3 cm	4,8 cm	RSZ-1503
645 cc	15,5 cm	15,9 cm	5,0 cm	RSZ-1553

CPG™ 332, Cohesive III™, Tall Height, Moderate *Plus* Projection

Vol.	Width	Height	Proj.	Catalog #
145 cc	9,0 cm	9,4 cm	3,8 cm	RSZ-0909
175 cc	9,5 cm	9,9 cm	4,0 cm	RSZ-0959
205 cc	10,0 cm	10,4 cm	4,2 cm	RSZ-1009
235 cc	10,5 cm	10,9 cm	4,4 cm	RSZ-1059
270 cc	11,0 cm	11,5 cm	4,7 cm	RSZ-1109
305 cc	11,5 cm	12,0 cm	4,9 cm	RSZ-1159
350 cc	12,0 cm	12,5 cm	5,1 cm	RSZ-1209
395 cc	12,5 cm	13,0 cm	5,3 cm	RSZ-1259
445 cc	13,0 cm	13,5 cm	5,5 cm	RSZ-1309
495 cc	13,5 cm	14,1 cm	5,7 cm	RSZ-1359
555 cc	14,0 cm	14,6 cm	5,9 cm	RSZ-1409
615 cc	14,5 cm	15,1 cm	6,1 cm	RSZ-1459
680 cc	15,0 cm	15,6 cm	6,3 cm	RSZ-1509

CPG™ 333, Cohesive III™, Tall Height, High Projection

Vol.	Width	Height	Proj.	Catalog #
180 cc	9,0 cm	9,4 cm	4,6 cm	RSZ-4090
215 cc	9,5 cm	9,9 cm	4,8 cm	RSZ-4095
250 cc	10,0 cm	10,4 cm	5,1 cm	RSZ-4100
290 cc	10,5 cm	10,9 cm	5,3 cm	RSZ-4105
330 cc	11,0 cm	11,5 cm	5,6 cm	RSZ-4110
380 cc	11,5 cm	12,0 cm	5,8 cm	RSZ-4115
430 cc	12,0 cm	12,5 cm	6,0 cm	RSZ-4120
485 cc	12,5 cm	13,0 cm	6,2 cm	RSZ-4125
545 cc	13,0 cm	13,5 cm	6,5 cm	RSZ-4130
610 cc	13,5 cm	14,1 cm	6,7 cm	RSZ-4135
680 cc	14,0 cm	14,6 cm	6,9 cm	RSZ-4140
755 cc	14,5 cm	15,1 cm	7,1 cm	RSZ-4145

Round Low Profile

Vol.	Diam.	Proj.	Catalog #
100 cc	9,4 cm	1,9 cm	RSZ-1100
125 cc	10,1 cm	2,0 cm	RSZ-1125
150 cc	10,7 cm	2,2 cm	RSZ-1150
175 cc	11,3 cm	2,3 cm	RSZ-1175
200 cc	11,8 cm	2,4 cm	RSZ-1200
225 cc	12,3 cm	2,5 cm	RSZ-1225
250 cc	12,7 cm	2,6 cm	RSZ-1250
275 cc	13,1 cm	2,7 cm	RSZ-1275
300 cc	13,5 cm	2,8 cm	RSZ-1300
325 cc	13,9 cm	2,9 cm	RSZ-1325
350 cc	14,2 cm	2,9 cm	RSZ-1350
375 cc	14,5 cm	2,9 cm	RSZ-1375
400 cc	14,8 cm	3,0 cm	RSZ-1400
450 cc	15,4 cm	3,1 cm	RSZ-1450
500 cc	16,0 cm	3,3 cm	RSZ-1500
550 cc	16,5 cm	3,4 cm	RSZ-1550
600 cc	17,0 cm	3,5 cm	RSZ-1600

Round Moderate Profile

Vol.	Diam.	Proj.	Catalog #
100 cc	9,3 cm	2,1 cm	RSZ-7100
125 cc	10,0 cm	2,2 cm	RSZ-7125
150 cc	10,3 cm	2,3 cm	RSZ-7150
175 cc	11,2 cm	2,4 cm	RSZ-7175
200 cc	11,7 cm	2,5 cm	RSZ-7200
225 cc	12,2 cm	2,6 cm	RSZ-7225
250 cc	12,3 cm	2,8 cm	RSZ-7250
275 cc	13,2 cm	2,9 cm	RSZ-7275
300 cc	13,5 cm	3,0 cm	RSZ-7300
325 cc	13,9 cm	3,0 cm	RSZ-7325
350 cc	14,2 cm	3,1 cm	RSZ-7350
375 cc	14,4 cm	3,2 cm	RSZ-7375
400 cc	14,5 cm	3,2 cm	RSZ-7400
450 cc	14,9 cm	3,4 cm	RSZ-7450
500 cc	15,2 cm	3,6 cm	RSZ-7500
550 cc	15,9 cm	3,8 cm	RSZ-7550
600 cc	16,5 cm	3,7 cm	RSZ-7600
700 cc	17,4 cm	3,9 cm	RSZ-7700
800 cc	18,2 cm	4,1 cm	RSZ-7800

CPG™ 321, Cohesive III™, Medium Height, Moderate Projection

Vol.	Width	Height	Proj.	Catalog #
120 cc	9,0 cm	8,5 cm	3,3 cm	RSZ-0908
135 cc	9,5 cm	8,9 cm	3,5 cm	RSZ-0958
155 cc	10,0 cm	9,4 cm	3,7 cm	RSZ-1008
180 cc	10,5 cm	9,9 cm	3,8 cm	RSZ-1058
215 cc	11,0 cm	10,3 cm	3,9 cm	RSZ-1108
245 cc	11,5 cm	10,8 cm	4,0 cm	RSZ-1158
280 cc	12,0 cm	11,3 cm	4,2 cm	RSZ-1208
315 cc	12,5 cm	11,8 cm	4,4 cm	RSZ-1258
355 cc	13,0 cm	12,2 cm	4,6 cm	RSZ-1308
395 cc	13,5 cm	12,7 cm	4,7 cm	RSZ-1358
440 cc	14,0 cm	13,2 cm	4,9 cm	RSZ-1408
480 cc	14,5 cm	13,6 cm	5,0 cm	RSZ-1458
530 cc	15,0 cm	14,1 cm	5,2 cm	RSZ-1508
640 cc	16,0 cm	15,0 cm	5,6 cm	RSZ-1608
775 cc	17,0 cm	16,0 cm	5,9 cm	RSZ-1708

CPG™ 322, Cohesive III™, Medium Height, Moderate *Plus* Projection

Vol.	Width	Height	Proj.	Catalog #
140 cc	9,0 cm	8,5 cm	3,8 cm	RSZ-0905
165 cc	9,5 cm	8,9 cm	4,0 cm	RSZ-0955
195 cc	10,0 cm	9,4 cm	4,2 cm	RSZ-1005
225 cc	10,5 cm	9,9 cm	4,4 cm	RSZ-1055
255 cc	11,0 cm	10,3 cm	4,7 cm	RSZ-1105
295 cc	11,5 cm	10,8 cm	4,9 cm	RSZ-1155
330 cc	12,0 cm	11,3 cm	5,1 cm	RSZ-1205
375 cc	12,5 cm	11,8 cm	5,3 cm	RSZ-1255
420 cc	13,0 cm	12,2 cm	5,5 cm	RSZ-1305
475 cc	13,5 cm	12,7 cm	5,7 cm	RSZ-1355
525 cc	14,0 cm	13,2 cm	5,9 cm	RSZ-1405
585 cc	14,5 cm	13,6 cm	6,1 cm	RSZ-1455
650 cc	15,0 cm	14,1 cm	6,3 cm	RSZ-1505

CPG™ 323, Cohesive III™, Medium Height, High Projection

Vol.	Width	Height	Proj.	Catalog #
185 cc	9,0 cm	8,5 cm	4,6 cm	RSZ-0902
195 cc	9,5 cm	8,9 cm	4,8 cm	RSZ-0952
225 cc	10,0 cm	9,4 cm	5,1 cm	RSZ-1002
260 cc	10,5 cm	9,9 cm	5,3 cm	RSZ-1052
300 cc	11,0 cm	10,3 cm	5,6 cm	RSZ-1102
345 cc	11,5 cm	10,8 cm	5,8 cm	RSZ-1152
390 cc	12,0 cm	11,3 cm	6,0 cm	RSZ-1202
440 cc	12,5 cm	11,8 cm	6,2 cm	RSZ-1252
495 cc	13,0 cm	12,2 cm	6,5 cm	RSZ-1302
555 cc	13,5 cm	12,7 cm	6,7 cm	RSZ-1352
620 cc	14,0 cm	13,2 cm	6,9 cm	RSZ-1402
685 cc	14,5 cm	13,6 cm	7,1 cm	RSZ-1452

Round Moderate *Plus* Profile

Vol.	Diam.	Proj.	Catalog #
100 cc	8,2 cm	2,7 cm	RSZ-1001
125 cc	8,9 cm	2,8 cm	RSZ-1251
150 cc	9,5 cm	2,9 cm	RSZ-1501
175 cc	10,0 cm	3,1 cm	RSZ-1751
200 cc	10,5 cm	3,2 cm	RSZ-2001
225 cc	10,9 cm	3,3 cm	RSZ-2251
250 cc	11,3 cm	3,4 cm	RSZ-2501
275 cc	11,7 cm	3,5 cm	RSZ-2751
300 cc	12,0 cm	3,6 cm	RSZ-3001
325 cc	12,3 cm	3,8 cm	RSZ-3251
350 cc	12,5 cm	3,9 cm	RSZ-3501
375 cc	12,8 cm	4,0 cm	RSZ-3751
400 cc	13,1 cm	4,0 cm	RSZ-4001
450 cc	13,6 cm	4,2 cm	RSZ-4501
500 cc	14,1 cm	4,3 cm	RSZ-5001
550 cc	14,6 cm	4,5 cm	RSZ-5501
600 cc	15,0 cm	4,6 cm	RSZ-6001
700 cc	15,8 cm	4,9 cm	RSZ-7001
800 cc	16,5 cm	5,1 cm	RSZ-8001

Round High Profile

Vol.	Diam.	Proj.	Catalog #
125 cc	8,3 cm	3,5 cm	RSZ-1254
150 cc	8,8 cm	3,7 cm	RSZ-1504
175 cc	9,3 cm	3,9 cm	RSZ-1754
200 cc	9,7 cm	4,0 cm	RSZ-2004
225 cc	10,1 cm	4,2 cm	RSZ-2254
250 cc	10,5 cm	4,3 cm	RSZ-2504
275 cc	10,8 cm	4,4 cm	RSZ-2754
300 cc	11,1 cm	4,5 cm	RSZ-3004
325 cc	11,4 cm	4,6 cm	RSZ-3254
350 cc	11,7 cm	4,8 cm	RSZ-3504
375 cc	12,0 cm	4,8 cm	RSZ-3754
400 cc	12,2 cm	5,0 cm	RSZ-4004
425 cc	12,5 cm	5,0 cm	RSZ-4254
450 cc	12,8 cm	5,1 cm	RSZ-4504
500 cc	13,2 cm	5,3 cm	RSZ-5004
550 cc	13,6 cm	5,5 cm	RSZ-5504
600 cc	14,0 cm	5,6 cm	RSZ-6004
650 cc	14,4 cm	5,7 cm	RSZ-6504
700 cc	14,8 cm	5,8 cm	RSZ-7004
800 cc	15,5 cm	6,0 cm	RSZ-8004

CPG™ 311^{NEW} Cohesive III™, Low Height, Moderate Projection

Vol.	Width	Height	Proj.	Catalog #
100 cc	9,0 cm	8,0 cm	3,0 cm	RSZ-0090
120 cc	9,5 cm	8,5 cm	3,2 cm	RSZ-0095
140 cc	10,0 cm	8,9 cm	3,4 cm	RSZ-0100
160 cc	10,5 cm	9,4 cm	3,6 cm	RSZ-0105
180 cc	11,0 cm	9,8 cm	3,7 cm	RSZ-0110
210 cc	11,5 cm	10,3 cm	3,9 cm	RSZ-0115
235 cc	12,0 cm	10,7 cm	4,1 cm	RSZ-0120
270 cc	12,5 cm	11,2 cm	4,2 cm	RSZ-0125
300 cc	13,0 cm	11,6 cm	4,4 cm	RSZ-0130
335 cc	13,5 cm	12,0 cm	4,6 cm	RSZ-0135
375 cc	14,0 cm	12,5 cm	4,7 cm	RSZ-0140
415 cc	14,5 cm	12,9 cm	4,9 cm	RSZ-0145
460 cc	15,0 cm	13,4 cm	5,1 cm	RSZ-0150
510 cc	15,5 cm	13,8 cm	5,2 cm	RSZ-0155
560 cc	16,0 cm	14,3 cm	5,4 cm	RSZ-0160
615 cc	16,5 cm	14,7 cm	5,6 cm	RSZ-0165

CPG™ 312, Cohesive III™, Low Height, Moderate *Plus* Projection

Vol.	Width	Height	Proj.	Catalog #
125 cc	9,0 cm	8,0 cm	3,8 cm	RSZ-0907
145 cc	9,5 cm	8,4 cm	4,0 cm	RSZ-0957
170 cc	10,0 cm	8,8 cm	4,2 cm	RSZ-1007
195 cc	10,5 cm	9,3 cm	4,4 cm	RSZ-1057
225 cc	11,0 cm	9,7 cm	4,7 cm	RSZ-1107
255 cc	11,5 cm	10,2 cm	4,9 cm	RSZ-1157
290 cc	12,0 cm	10,6 cm	5,1 cm	RSZ-1207
330 cc	12,5 cm	11,1 cm	5,3 cm	RSZ-1257
370 cc	13,0 cm	11,5 cm	5,5 cm	RSZ-1307
415 cc	13,5 cm	12,0 cm	5,7 cm	RSZ-1357
465 cc	14,0 cm	12,4 cm	5,9 cm	RSZ-1407
515 cc	14,5 cm	12,8 cm	6,1 cm	RSZ-1457
570 cc	15,0 cm	13,3 cm	6,3 cm	RSZ-1507
690 cc	16,0 cm	14,2 cm	6,8 cm	RSZ-1607

CPG™ 313^{NEW} Cohesive III™, Low Height, High Projection

Vol.	Width	Height	Proj.	Catalog #
130 cc	9,0 cm	8,1 cm	4,4 cm	RSZ-0906
155 cc	9,5 cm	8,6 cm	4,5 cm	RSZ-0956
180 cc	10,0 cm	9,0 cm	4,6 cm	RSZ-1006
210 cc	10,5 cm	9,5 cm	4,8 cm	RSZ-1056
240 cc	11,0 cm	9,9 cm	4,9 cm	