

Smooth Round Saline Implant Family



Style 2000 Moderate *Plus* Profile



Style 3000 High Profile



Style 1600 Moderate Profile

- Available in Moderate, Moderate *Plus* and High Profile

- Smooth Surface Shell

- Anterior Diaphragm Valve

- Intrasurgical Size Adjustability (25-75 cc for Moderate, 35-160 cc for Moderate Plus and 30-120 cc for High Profile)

Please ensure that all patients receive a copy of the Informed Decision brochure prior to deciding on surgery for saline breast implants. For the latest information on saline breast implants, please refer to the Product Insert provided with each implant or on our website www.mentorcorp.com.

Smooth Round Saline Breast Implant Family

Smooth Round Moderate Profile

NOMINAL FILL			MAXIMUM FILL			Catalog #	Sterile Sizer
Vol.	Diam.	Proj.	Vol.	Diam.	Proj.		
125 cc	9.5 cm	3.0 cm	150 cc	9.1 cm	3.4 cm	350-1610	351-125SZ
150 cc	10.0 cm	3.1 cm	175 cc	9.6 cm	3.7 cm	350-1615	
175 cc	10.6 cm	3.3 cm	200 cc	10.0 cm	3.9 cm	350-1620	351-175SZ
200 cc	11.0 cm	3.4 cm	225 cc	10.4 cm	3.9 cm	350-1625	
225 cc	11.5 cm	3.5 cm	250 cc	10.9 cm	4.0 cm	350-1630	351-225SZ
250 cc	11.9 cm	3.6 cm	275 cc	11.2 cm	4.1 cm	350-1635	
275 cc	12.3 cm	3.7 cm	300 cc	11.7 cm	4.2 cm	350-1640	351-275SZ
300 cc	12.6 cm	3.7 cm	325 cc	12.1 cm	4.3 cm	350-1645	
325 cc	13.0 cm	3.8 cm	375 cc	12.3 cm	4.5 cm	350-1650	351-325SZ
350 cc	13.3 cm	3.9 cm	400 cc	12.6 cm	4.7 cm	350-1655	
375 cc	13.6 cm	4.0 cm	425 cc	12.8 cm	4.8 cm	350-1660	351-375SZ
425 cc	14.2 cm	4.1 cm	475 cc	13.5 cm	4.9 cm	350-1670	351-425SZ
475 cc	14.8 cm	4.2 cm	525 cc	14.2 cm	4.9 cm	350-1680	351-475SZ
525 cc	15.0 cm	4.2 cm	575 cc	14.5 cm	5.1 cm	350-1685	
575 cc	15.0 cm	4.5 cm	625 cc	14.8 cm	5.1 cm	350-1690	
625 cc	15.2 cm	4.6 cm	700 cc	14.9 cm	5.6 cm	350-1695	
700 cc	15.6 cm	4.9 cm	775 cc	14.9 cm	5.8 cm	350-1697	

Smooth Round Moderate *Plus* Profile

NOMINAL FILL			MID-RANGE FILL		MAXIMUM FILL		Catalog #	Sterile Sizer
Vol.	Diam.	Proj.	Diam.	Proj.	Diam.	Proj.		
175+35 cc	9.5 cm	3.5 cm	9.4 cm	3.9 cm	9.3 cm	4.2 cm	350-2175	351-2175SZ
200+40 cc	10.0 cm	3.7 cm	9.9 cm	4.0 cm	9.7 cm	4.4 cm	350-2200	351-2200SZ
225+45 cc	10.4 cm	3.8 cm	10.3 cm	4.2 cm	10.1 cm	4.6 cm	350-2225	351-2225SZ
250+50 cc	10.8 cm	4.0 cm	10.7 cm	4.3 cm	10.5 cm	4.7 cm	350-2250	351-2250SZ
275+55 cc	11.0 cm	4.1 cm	10.9 cm	4.5 cm	10.8 cm	4.9 cm	350-2275	351-2275SZ
300+60 cc	11.5 cm	4.3 cm	11.4 cm	4.6 cm	11.2 cm	5.0 cm	350-2300	351-2300SZ
325+65 cc	11.9 cm	4.4 cm	11.7 cm	4.7 cm	11.5 cm	5.2 cm	350-2325	351-2325SZ
350+70 cc	12.1 cm	4.5 cm	12.0 cm	4.8 cm	11.7 cm	5.3 cm	350-2350	351-2350SZ
375+75 cc	12.3 cm	4.6 cm	12.2 cm	4.9 cm	12.0 cm	5.4 cm	350-2375	351-2375SZ
400+80 cc	12.6 cm	4.7 cm	12.5 cm	5.0 cm	12.3 cm	5.5 cm	350-2400	351-2400SZ
425+85 cc	12.9 cm	4.8 cm	12.8 cm	5.1 cm	12.5 cm	5.6 cm	350-2425	351-2425SZ
450+90 cc	13.0 cm	4.9 cm	12.9 cm	5.2 cm	12.8 cm	5.7 cm	350-2450	351-2450SZ
475+95 cc	13.3 cm	5.0 cm	13.2 cm	5.3 cm	13.0 cm	5.8 cm	350-2475	351-2475SZ
500+100 cc	13.6 cm	5.1 cm	13.5 cm	5.4 cm	13.2 cm	5.9 cm	350-2500	351-2500SZ
550+110 cc	14.0 cm	5.3 cm	13.9 cm	5.6 cm	13.7 cm	6.1 cm	350-2550	351-2550SZ
600+120 cc	14.5 cm	5.5 cm	14.3 cm	5.8 cm	14.1 cm	6.3 cm	350-2600	351-2600SZ
650+130 cc	14.8 cm	5.6 cm	14.7 cm	5.9 cm	14.4 cm	6.5 cm	350-2650	351-2650SZ
700+140 cc	15.2 cm	5.8 cm	15.1 cm	6.1 cm	14.8 cm	6.6 cm	350-2700	351-2700SZ
750+150 cc	15.6 cm	5.9 cm	15.4 cm	6.2 cm	15.1 cm	6.8 cm	350-2750	351-2750SZ
800+160 cc	16.1 cm	6.0 cm	15.8 cm	6.3 cm	15.5 cm	7.0 cm	350-2800	351-2800SZ

Smooth Round High Profile

Vol.	NOMINAL FILL		MID-RANGE FILL		MAXIMUM FILL		Catalog #	Sterile Sizer
	Diam.	Proj.	Diam.	Proj.	Diam.	Proj.		
170+30 cc	9.0 cm	4.0 cm	8.8 cm	4.5 cm	8.8 cm	4.9 cm	350-3170	
190+35 cc	9.3 cm	4.1 cm	9.2 cm	4.6 cm	9.1 cm	5.1 cm	350-3190	351-3190SZ
210+40 cc	9.6 cm	4.2 cm	9.5 cm	4.8 cm	9.5 cm	5.3 cm	350-3210	
230+45 cc	10.0 cm	4.3 cm	9.8 cm	4.9 cm	9.8 cm	5.5 cm	350-3230	351-3230SZ
250+50 cc	10.2 cm	4.5 cm	10.1 cm	4.9 cm	10.0 cm	5.6 cm	350-3250	
270+55 cc	10.4 cm	4.6 cm	10.5 cm	5.0 cm	10.2 cm	5.8 cm	350-3270	351-3270SZ
290+60 cc	10.8 cm	4.7 cm	10.7 cm	5.1 cm	10.5 cm	5.9 cm	350-3290	351-3290SZ
310+65 cc	11.0 cm	4.8 cm	10.9 cm	5.3 cm	10.7 cm	6.0 cm	350-3310	351-3310SZ
330+70 cc	11.3 cm	4.8 cm	11.2 cm	5.4 cm	11.0 cm	6.2 cm	350-3330	351-3330SZ
380+70 cc	11.7 cm	5.2 cm	11.6 cm	5.7 cm	11.4 cm	6.4 cm	350-3380	351-3380SZ
420+80 cc	12.0 cm	5.4 cm	12.0 cm	5.8 cm	11.7 cm	6.7 cm	350-3420	351-3420SZ
460+90 cc	12.4 cm	5.5 cm	12.6 cm	6.0 cm	12.1 cm	6.9 cm	350-3460	
500+100 cc	12.8 cm	5.6 cm	12.9 cm	6.1 cm	12.4 cm	7.1 cm	350-3500	351-3500SZ
560+115 cc	13.4 cm	5.7 cm	13.3 cm	6.4 cm	13.1 cm	7.3 cm	350-3560	
630+120 cc	13.8 cm	5.9 cm	13.7 cm	6.7 cm	13.4 cm	7.4 cm	350-3630	351-3630SZ

Diameter Comparison



Note: Individual implant dimensions may vary slightly in products of this type. Not all units will conform exactly to the dimensions noted above.

International Operations

Mentor Medical Systems, Australia, Pty. Ltd. Tel: +61 2 9251 8111 Toll Free: +1 800 338 160 Fax: +61 2 9251 8333

Mentor Medical Systems, Canada, Inc. Tel: +1 905 725 7763 Toll Free: +1 800 668 6069 Fax: +1 905 725 7340

Mentor Deutschland GmbH Tel: +49 811 600 50 0 Fax: +49 811 600 50 11

Mentor Medical Systems France, S.A. Tel: +33 1 30 07 07 60 Fax: +33 1 30 07 07 67

Mentor Medical Systems Iberica, S.L. Tel: +34 91 562 2700 Fax: +34 91 728 3367

Mentor Medical Italia, S.r.l. Tel: +39 0 2 880 7761 Fax: +39 0 2 7209 3710

Mentor Medical Systems Ltd., U.K. Tel: +44 1635 511800 Fax: +44 1635 511810



201 Mentor Drive
Santa Barbara, CA 93111 USA

(805) 879-6000

Customer Service:
Tel: (800) 235-5731
Fax: (805) 967-7108

www.mentorcorp.com

Manufacturer:
Mentor
3041 Skyway Circle North
Irving, TX 75038 USA

European Representative:
Mentor Medical Systems B.V.
2333 CL Leiden, The Netherlands