

Smooth Tissue Expander Family




- Available in Round, Rectangular, Crescent, and Elliptical Shapes
- Dacron-reinforced Base Provides a Firm Platform for Directional Expansion
- Bottom-loading Fill Design Minimizes Potential Side-wall Leakage
- Remote Injection Dome

Not for distribution in the U.S.

Smooth Tissue Expander Family


Round Tissue Expander

Volume	Diameter	Projection	Catalog Number
400 cc	11,3 cm	5,9 cm	350-4305M
550 cc	13,4 cm	5,8 cm	350-4307M
700 cc	14,6 cm	6,7 cm	350-4304M
850 cc	15,0 cm	7,2 cm	350-4308M
1000 cc	15,5 cm	7,8 cm	350-4306M




Rectangle Tissue Expander

Volume	Length	Width	Projection	Catalog Number
50 cc	6,9 cm	3,2 cm	3,2 cm	350-4309M
100 cc	6,6 cm	5,1 cm	5,6 cm	350-4301M
250 cc	10,4 cm	7,3 cm	5,6 cm	350-4302M
400 cc	10,6 cm	9,3 cm	6,6 cm	350-4303M
550 cc	13,7 cm	7,1 cm	7,7 cm	350-4310M
700 cc	15,0 cm	8,0 cm	7,9 cm	350-4311M



Crescent Tissue Expander


Volume	Length	Width	Projection	Catalog Number
50 cc	9,0 cm	5,7 cm	2,7 cm	350-4312M
100 cc	11,0 cm	6,9 cm	3,7 cm	350-4313M
200 cc	14,6 cm	7,5 cm	4,3 cm	350-4314M
300 cc	16,5 cm	10,3 cm	4,7 cm	350-4315M
400 cc	18,2 cm	11,4 cm	5,2 cm	350-4316M
500 cc	19,6 cm	12,4 cm	5,5 cm	350-4317M



Elliptical Tissue Expander

TEMPORARY OVEREXPANSION

Device Volume	Length	Width	Normal Projection	Projection	Max. Volume	Catalog Number
25 cc	6,0 cm	4,0 cm	2,3 cm	4,5 cm	75 cc	350-5301M
50 cc	8,0 cm	5,0 cm	2,5 cm	5,3 cm	150 cc	350-5302M
75 cc	10,0 cm	6,0 cm	2,1 cm	6,0 cm	225 cc	350-5303M
150 cc	12,0 cm	7,0 cm	3,2 cm	7,2 cm	450 cc	350-5304M
250 cc	14,0 cm	8,0 cm	4,5 cm	9,0 cm	750 cc	350-5305M
125 cc	15,0 cm	5,0 cm	3,3 cm	6,5 cm	375 cc	350-5306M



Accessory Items

Description	Catalog Number
Injection Dome Pack (contains micro and standard domes)	350-DOMPK
21-Gauge Winged Infusion Set (each)	200558-001
Tissue Expander Templates	LIT308

Tissue Expanders are packaged with:

- Standard remote injection dome
- Stainless steel connector (to allow customization of tube length)
- 23-gauge butterfly needle

PLEASE NOTE: Individual implant dimensions may vary slightly in products of this type. Not all units conform exactly to the dimensions noted above.

For detailed indications, contraindications, warnings and precautions associated with the use of this device, please refer to the Product Insert Data Sheet.

Important Safety Information:

Contour Profile® Tissue Expanders and Smooth Round Tissue Expanders with Remote Dome are used for breast reconstruction following mastectomy. These expanders are intended for temporary subcutaneous or submuscular implantation and are not intended for use beyond six months. Do not use the Contour Profile® Tissue Expander and the Smooth Round Tissue Expander with Remote Dome in patients where an MRI may be needed. The device could be moved by the MRI causing pain or displacement, potentially resulting in a revision surgery. The incidence of extrusion of the expander has been shown to increase when the expander has been placed in injured areas.

For detailed indications, contraindications, warning and precautions associated with the use of Tissue Expanders, please refer to the Product Insert Data Sheet provided with each product, or review the Important Safety Information provided at www.mentorcorp.com.



201 Mentor Drive
Santa Barbara, CA 93111 USA
Tel: +1 805 879 6000
Fax: +1 805 967 7108
www.mentorcorp.com

Manufacturer:
Mentor
3041 Skyway Circle North
Irving, TX 75038 USA

European Representative:
Mentor Medical Systems B.V.
2333 CL Leiden, The Netherlands

International Operations

Mentor Medical Systems, Australia, Pty. Ltd. Tel: +61 2 9251 8111 Toll Free: +1 800 338 160 Fax: +61 2 9251 8333
Mentor Medical Systems, Canada, Inc. Tel: +1 905 725 7763 Toll Free: +1 800 668 6069 Fax: +1 905 725 7340
Mentor Deutschland GmbH Tel: +49 811 600 50 0 Fax: +49 811 600 50 11
Mentor Medical Systems Iberica, S.L. Tel: +34 91 562 2700 Fax: +34 91 728 3367
Mentor Medical Italia, S.r.l. Tel: +39 0 2 880 7761 Fax: +39 0 2 7209 3710
Mentor Medical Systems Ltd., U.K. Tel: +44 1635 511800 Fax: +44 1635 511810
Perouse Plastie, S.A.S., France Tel: +33 3 44 08 26 26 Fax: +33 3 44 08 26 01