


Siltex™ Round Gel Implant Family

Cohesive II™




Moderate *Plus* Profile



High Profile



Moderate Profile



Cohesive II™ gel is a slightly firmer gel, for surgeons wanting a firmer feeling implant. This gel is used in Round Moderate Profile, Moderate Plus Profile, and High Profile.

- Siltex™ Textured Surface, available in Moderate, Moderate Plus, and High Profile

- Offered in a Wide Range of Sizes

- Features Mentor's Proprietary Cohesive Gel Formulation

- Shell Barrier Layer Designed to Reduce Gel Bleed

- Lifetime Product Replacement Policy

Not for distribution in the U.S.

Siltex™ Round Gel Breast Implant Family

Siltex™ Round Moderate Profile, Cohesive II™

Vol. cc	Diam. cm	Proj. cm	Catalog Number	Resterilizable Gel Sizer*
100	9,1	2,6	354-1009	RSZ-7100
125	9,3	2,8	354-1259	RSZ-7125
150	10,1	2,4	354-1509	RSZ-7150
175	10,5	2,9	354-1759	RSZ-7175
200	11,2	2,8	354-2009	RSZ-7200
225	11,5	3,0	354-2259	RSZ-7225
250	11,6	3,5	354-2509	RSZ-7250
275	12,6	3,1	354-2759	RSZ-7275
300	12,6	3,6	354-3009	RSZ-7300
325	13,1	3,6	354-3259	RSZ-7325
350	13,4	3,5	354-3509	RSZ-7350
375	13,4	3,9	354-3759	RSZ-7375
400	13,7	3,8	354-4009	RSZ-7400
450	13,7	4,6	354-4509	RSZ-7450
500	14,2	4,8	354-5009	RSZ-7500
550	14,5	4,5	354-5509	RSZ-7550
600	15,3	4,6	354-6009	RSZ-7600
700	16,3	5,0	354-7009	RSZ-7700
800	16,7	5,4	354-8009	RSZ-7800

offered in 19 sizes

Siltex™ Round Moderate *Plus* Profile, Cohesive II™

Vol. cc	Diam. cm	Proj. cm	Catalog Number	Resterilizable Gel Sizer*
100	8,1	2,7	324-5100	RSZ-1001
125	8,8	2,9	324-5125	RSZ-1251
150	9,4	3,0	324-5150	RSZ-1501
175	10,0	3,2	324-5175	RSZ-1751
200	10,5	3,3	324-5200	RSZ-2001
225	10,9	3,5	324-5225	RSZ-2251
250	11,3	3,6	324-5250	RSZ-2501
275	11,7	3,7	324-5275	RSZ-2751
300	12,0	3,7	324-5300	RSZ-3001
325	12,3	3,8	324-5325	RSZ-3251
350	12,6	3,8	324-5350	RSZ-3501
375	12,9	3,9	324-5375	RSZ-3751
400	13,2	4,0	324-5400	RSZ-4001
450	13,7	4,1	324-5450	RSZ-4501
500	14,1	4,2	324-5500	RSZ-5001
550	14,4	4,4	324-5550	RSZ-5501
600	14,7	4,5	324-5600	RSZ-6001
700	15,7	4,8	324-5700	RSZ-7001
800	16,6	5,0	324-5800	RSZ-8001

offered in 19 sizes

Siltex™ Round High Profile, Cohesive II™


Vol. cc	Diam. cm	Proj. cm	Catalog Number	Resterilizable Gel Sizer*
125	8,4	3,6	324-4125	RSZ-1254
150	8,9	3,8	324-4150	RSZ-1504
175	9,4	4,0	324-4175	RSZ-1754
200	9,9	4,1	324-4200	RSZ-2004
225	10,2	4,3	324-4225	RSZ-2254
250	10,5	4,5	324-4250	RSZ-2504
275	10,9	4,6	324-4275	RSZ-2754
300	11,1	4,7	324-4300	RSZ-3004
325	11,5	4,8	324-4325	RSZ-3254
350	11,7	4,9	324-4350	RSZ-3504
375	12,0	5,0	324-4375	RSZ-3754
400	12,3	5,1	324-4400	RSZ-4004
425	12,5	5,2	324-4425	RSZ-4254
450	12,7	5,2	324-4450	RSZ-4504
500	13,2	5,4	324-4500	RSZ-5004
550	13,5	5,6	324-4550	RSZ-5504
600	14,0	5,7	324-4600	RSZ-6004
650	14,3	5,8	324-4650	RSZ-6504
700	14,7	6,0	324-4700	RSZ-7004
800	15,4	6,3	324-4800	RSZ-8004

offered in 20 sizes

Note: Individual implant dimensions may vary slightly in products of this type. Not all units will conform exactly to the dimensions noted above.

* Resterilizable up to 10 times.

225 cc Comparison



INTERNATIONAL OPERATIONS

Mentor Medical Systems, Australia, Pty. Ltd. Tel: +61 2 9251 8111 Toll Free: +1 800 338 160 Fax: +61 2 9251 8333

Mentor Medical Systems, Canada, Inc. Tel: +1 905 725 7763 Toll Free: +1 800 668 6069 Fax: +1 905 725 7340

Mentor Deutschland GmbH Tel: +49 811 600 50 0 Fax: +49 811 600 50 11

Mentor Medical Systems Iberica, S.L. Tel: +34 91 562 2700 Fax: +34 91 728 3367

Mentor Medical Italia, S.r.l. Tel: +39 0 2 880 7761 Fax: +39 0 2 7209 3710

Mentor Medical Systems Ltd., U.K. Tel: +44 1635 511800 Fax: +44 1635 511810

Perouse Plastique, S.A.S. Tel: +33 3 44 08 26 26 Fax: +33 3 44 08 26 01



201 Mentor Drive
Santa Barbara, CA 93111 USA

+1 (805) 879-6000
Fax: +1 (805) 967-7108
www.mentorcorp.com

Manufacturer:
Mentor Medical Systems B.V.
Zernikedreef 2
2333 CL Leiden
The Netherlands