


Smooth Round Gel Implant Family

Cohesive I™



Low Profile



Moderate Profile



Moderate *Plus* Profile



Ultra High Profile



High Profile



Cohesive I™ gel is the standard gel used in Mentor's implants. This is Mentor's softest gel and has been preferred for years in Round Low, Moderate, Moderate *Plus*, High Profile, Ultra High Profile, and Becker implants.

- Smooth Surface, available in Low, Moderate, Moderate Plus, High, and Ultra High Profile
- Offered in a Wide Range of Sizes
- Features Mentor's Proprietary Cohesive Gel Formulation
- Shell Barrier Layer Designed to Reduce Gel Bleed
- Lifetime Product Replacement Policy

Not for distribution in the U.S.



Smooth Round Gel Breast Implant Family

Smooth Round Low Profile, Cohesive I™

| Vol. | Diam. | Proj. | Catalog # | Restorable Gel Sizer* |
|--------|---------|--------|------------|--------------------------|
| 100 cc | 9,4 cm | 1,9 cm | 310JSP1100 | RSZ-1100 |
| 125 cc | 10,1 cm | 2,0 cm | 310JSP1125 | RSZ-1125 |
| 150 cc | 10,7 cm | 2,2 cm | 310JSP1150 | RSZ-1150 |
| 175 cc | 11,3 cm | 2,3 cm | 310JSP1175 | RSZ-1175 |
| 200 cc | 11,8 cm | 2,4 cm | 310JSP1200 | RSZ-1200 |
| 225 cc | 12,3 cm | 2,5 cm | 310JSP1225 | RSZ-1225 |
| 250 cc | 12,7 cm | 2,6 cm | 310JSP1250 | RSZ-1250 |
| 275 cc | 13,1 cm | 2,7 cm | 310JSP1275 | RSZ-1275 |
| 300 cc | 13,5 cm | 2,7 cm | 310JSP1300 | RSZ-1300 |
| 325 cc | 13,9 cm | 2,8 cm | 310JSP1325 | RSZ-1325 |
| 350 cc | 14,2 cm | 2,9 cm | 310JSP1350 | RSZ-1350 |
| 375 cc | 14,5 cm | 2,9 cm | 310JSP1375 | RSZ-1375 |
| 400 cc | 14,8 cm | 3,0 cm | 310JSP1400 | RSZ-1400 |
| 450 cc | 15,4 cm | 3,1 cm | 310JSP1450 | RSZ-1450 |
| 500 cc | 16,0 cm | 3,3 cm | 310JSP1500 | RSZ-1500 |
| 550 cc | 16,5 cm | 3,4 cm | 310JSP1550 | RSZ-1550 |
| 600 cc | 17,0 cm | 3,5 cm | 310JSP1600 | RSZ-1600 |

offered in 17 sizes

Smooth Round Moderate Profile, Cohesive I™

| Vol. | Diam. | Proj. | Catalog # | Restorable Gel Sizer* |
|--------|---------|--------|-------------|--------------------------|
| 100 cc | 9,3 cm | 2,1 cm | 350-7100 BC | RSZ-7100 |
| 125 cc | 10,0 cm | 2,2 cm | 350-7125 BC | RSZ-7125 |
| 150 cc | 10,3 cm | 2,3 cm | 350-7150 BC | RSZ-7150 |
| 175 cc | 11,2 cm | 2,4 cm | 350-7175 BC | RSZ-7175 |
| 200 cc | 11,7 cm | 2,5 cm | 350-7200 BC | RSZ-7200 |
| 225 cc | 12,2 cm | 2,6 cm | 350-7225 BC | RSZ-7225 |
| 250 cc | 12,3 cm | 2,8 cm | 350-7250 BC | RSZ-7250 |
| 275 cc | 13,2 cm | 2,9 cm | 350-7275 BC | RSZ-7275 |
| 300 cc | 13,5 cm | 3,0 cm | 350-7300 BC | RSZ-7300 |
| 325 cc | 13,9 cm | 3,0 cm | 350-7325 BC | RSZ-7325 |
| 350 cc | 14,2 cm | 3,1 cm | 350-7350 BC | RSZ-7350 |
| 375 cc | 14,4 cm | 3,2 cm | 350-7375 BC | RSZ-7375 |
| 400 cc | 14,5 cm | 3,2 cm | 350-7400 BC | RSZ-7400 |
| 450 cc | 14,9 cm | 3,4 cm | 350-7450 BC | RSZ-7450 |
| 500 cc | 15,2 cm | 3,6 cm | 350-7500 BC | RSZ-7500 |
| 550 cc | 15,9 cm | 3,6 cm | 350-7550 BC | RSZ-7550 |
| 600 cc | 16,5 cm | 3,7 cm | 350-7600 BC | RSZ-7600 |
| 700 cc | 17,4 cm | 3,9 cm | 350-7700 BC | RSZ-7700 |
| 800 cc | 18,2 cm | 4,1 cm | 350-7800 BC | RSZ-7800 |


offered in 19 sizes

Smooth Round Moderate Plus Profile, Cohesive I™

| Vol. | Diam. | Proj. | Catalog # | Restorable Gel Sizer* |
|--------|---------|--------|-------------|--------------------------|
| 100 cc | 8,2 cm | 2,7 cm | 350-1001 BC | RSZ-1001 |
| 125 cc | 8,9 cm | 2,8 cm | 350-1251 BC | RSZ-1251 |
| 150 cc | 9,5 cm | 2,9 cm | 350-1501 BC | RSZ-1501 |
| 175 cc | 10,0 cm | 3,1 cm | 350-1751 BC | RSZ-1751 |
| 200 cc | 10,5 cm | 3,2 cm | 350-2001 BC | RSZ-2001 |
| 225 cc | 10,9 cm | 3,3 cm | 350-2251 BC | RSZ-2251 |
| 250 cc | 11,3 cm | 3,4 cm | 350-2501 BC | RSZ-2501 |
| 275 cc | 11,7 cm | 3,5 cm | 350-2751 BC | RSZ-2751 |
| 300 cc | 12,0 cm | 3,6 cm | 350-3001 BC | RSZ-3001 |
| 325 cc | 12,3 cm | 3,8 cm | 350-3251 BC | RSZ-3251 |
| 350 cc | 12,5 cm | 3,9 cm | 350-3501 BC | RSZ-3501 |
| 375 cc | 12,8 cm | 4,0 cm | 350-3751 BC | RSZ-3751 |
| 400 cc | 13,1 cm | 4,0 cm | 350-4001 BC | RSZ-4001 |
| 450 cc | 13,6 cm | 4,2 cm | 350-4501 BC | RSZ-4501 |
| 500 cc | 14,1 cm | 4,3 cm | 350-5001 BC | RSZ-5001 |
| 550 cc | 14,6 cm | 4,5 cm | 350-5501 BC | RSZ-5501 |
| 600 cc | 15,0 cm | 4,6 cm | 350-6001 BC | RSZ-6001 |
| 700 cc | 15,8 cm | 4,9 cm | 350-7001 BC | RSZ-7001 |
| 800 cc | 16,5 cm | 5,1 cm | 350-8001 BC | RSZ-8001 |

offered in 19 sizes

350 cc Comparison



Note: Individual implant dimensions may vary slightly in products of this type. Not all units will conform exactly to the dimensions noted above.

* Resterilizable up to 10 times.

Smooth Round High Profile, Cohesive I™

| Vol. | Diam. | Proj. | Catalog # | Restorable Gel Sizer* |
|--------|---------|--------|-------------|--------------------------|
| 125 cc | 8,3 cm | 3,5 cm | 350-1254 BC | RSZ-1254 |
| 150 cc | 8,8 cm | 3,7 cm | 350-1504 BC | RSZ-1504 |
| 175 cc | 9,3 cm | 3,9 cm | 350-1754 BC | RSZ-1754 |
| 200 cc | 9,7 cm | 4,0 cm | 350-2004 BC | RSZ-2004 |
| 225 cc | 10,1 cm | 4,2 cm | 350-2254 BC | RSZ-2254 |
| 250 cc | 10,5 cm | 4,3 cm | 350-2504 BC | RSZ-2504 |
| 275 cc | 10,8 cm | 4,4 cm | 350-2754 BC | RSZ-2754 |
| 300 cc | 11,1 cm | 4,5 cm | 350-3004 BC | RSZ-3004 |
| 325 cc | 11,4 cm | 4,6 cm | 350-3254 BC | RSZ-3254 |
| 350 cc | 11,7 cm | 4,8 cm | 350-3504 BC | RSZ-3504 |
| 375 cc | 12,0 cm | 4,8 cm | 350-3754 BC | RSZ-3754 |
| 400 cc | 12,2 cm | 5,0 cm | 350-4004 BC | RSZ-4004 |
| 425 cc | 12,5 cm | 5,0 cm | 350-4254 BC | RSZ-4254 |
| 450 cc | 12,8 cm | 5,1 cm | 350-4504 BC | RSZ-4504 |
| 500 cc | 13,2 cm | 5,3 cm | 350-5004 BC | RSZ-5004 |
| 550 cc | 13,6 cm | 5,5 cm | 350-5504 BC | RSZ-5504 |
| 600 cc | 14,0 cm | 5,6 cm | 350-6004 BC | RSZ-6004 |
| 650 cc | 14,4 cm | 5,7 cm | 350-6504 BC | RSZ-6504 |
| 700 cc | 14,8 cm | 5,8 cm | 350-7004 BC | RSZ-7004 |
| 800 cc | 15,5 cm | 6,0 cm | 350-8004 BC | RSZ-8004 |

offered in 20 sizes

Smooth Round Ultra High Profile, Cohesive I™

| Vol. | Diam. | Proj. | Catalog # | Restorable Gel Sizer* |
|--------|---------|--------|-------------|--------------------------|
| 135 cc | 7,8 cm | 4,1 cm | 350-5135 BC | RSZ-5135 |
| 160 cc | 8,2 cm | 4,3 cm | 350-5160 BC | RSZ-5160 |
| 185 cc | 8,4 cm | 4,4 cm | 350-5185 BC | RSZ-5185 |
| 215 cc | 8,7 cm | 4,5 cm | 350-5215 BC | RSZ-5215 |
| 240 cc | 9,0 cm | 4,7 cm | 350-5240 BC | RSZ-5240 |
| 270 cc | 9,3 cm | 4,8 cm | 350-5270 BC | RSZ-5270 |
| 295 cc | 9,7 cm | 5,0 cm | 350-5295 BC | RSZ-5295 |
| 320 cc | 9,8 cm | 5,1 cm | 350-5320 BC | RSZ-5320 |
| 350 cc | 10,1 cm | 5,2 cm | 350-5350 BC | RSZ-5350 |
| 375 cc | 10,4 cm | 5,3 cm | 350-5375 BC | RSZ-5375 |
| 400 cc | 10,6 cm | 5,4 cm | 350-5400 BC | RSZ-5400 |
| 430 cc | 10,9 cm | 5,6 cm | 350-5430 BC | RSZ-5430 |
| 455 cc | 11,2 cm | 5,7 cm | 350-5455 BC | RSZ-5455 |
| 480 cc | 11,4 cm | 5,8 cm | 350-5480 BC | RSZ-5480 |
| 535 cc | 12,0 cm | 6,1 cm | 350-5535 BC | RSZ-5535 |
| 590 cc | 12,5 cm | 6,3 cm | 350-5590 BC | RSZ-5590 |

offered in 16 sizes



201 Mentor Drive
Santa Barbara, CA 93111 USA
Tel: +1 (805) 879-6000
Fax: +1 (805) 967-7108
www.mentorcorp.com

Manufacturer:
Mentor Medical Systems B.V.
Zernikedreef 2
2333 CL Leiden
The Netherlands

INTERNATIONAL OPERATIONS

Mentor Medical Systems, Australia, Pty. Ltd. Tel: +61 2 9251 8111 Toll Free: +1 800 338 160 Fax: +61 2 9251 8333

Mentor Medical Systems, Canada, Inc. Tel: +1 905 725 7763 Toll Free: +1 800 668 6069 Fax: +1 905 725 7340

Mentor Deutschland GmbH Tel: +49 811 600 50 0 Fax: +49 811 600 50 11

Mentor Medical Systems Iberica, S.L. Tel: +34 91 562 2700 Fax: +34 91 728 3367

Mentor Medical Italia, S.r.l. Tel: +39 0 2 880 7761 Fax: +39 0 2 7209 3710

Mentor Medical Systems Ltd., U.K. Tel: +44 1635 511800 Fax: +44 1635 511810

Perouse Plastie, S.A.S. Tel: +33 3 44 08 26 26 Fax: +33 3 44 08 26 01