

Spectrum® Saline Implant Family



- Posterior Spectrum Valve

- Preplaced Fill-tube

- Remote Injection Dome for Optional Postoperative Adjustability

- Available in Three Styles
 - Round with Smooth Surface
 - Round with Siltex Textured Surface
 - Contour Profile with Siltex Textured Surface

Please ensure that all patients receive a copy of the Informed Decision brochure prior to deciding on surgery for saline breast implants. Please refer to new Product Insert Data Sheet provided with each implant for the latest information on saline breast implants.

Spectrum® Saline Breast Implant Family

Smooth Round Spectrum®

MINIMUM FILL

MAXIMUM FILL

Device Volume	Catalog Number	Temp. Min. Volume	Final Min. Volume	Min. Diameter	Min. Projection	Final Max. Volume	Max. Diameter	Max. Projection
125 cc	350-1410	105 cc	125 cc	9.4 cm	2.7 cm	150 cc	9.3 cm	3.2 cm
175 cc	350-1420	150 cc	175 cc	10.3 cm	3.1 cm	210 cc	10.2 cm	3.8 cm
225 cc	350-1430	190 cc	225 cc	10.8 cm	3.5 cm	270 cc	10.9 cm	4.3 cm
275 cc	350-1440	230 cc	275 cc	11.6 cm	3.8 cm	330 cc	11.5 cm	4.6 cm
325 cc	350-1450	275 cc	325 cc	12.4 cm	4.0 cm	390 cc	12.1 cm	4.9 cm
375 cc	350-1460	320 cc	375 cc	13.5 cm	4.1 cm	450 cc	13.0 cm	4.9 cm
425 cc	350-1470	360 cc	425 cc	14.0 cm	4.2 cm	510 cc	13.7 cm	4.8 cm
475 cc	350-1480	405 cc	475 cc	14.3 cm	4.4 cm	570 cc	13.8 cm	5.4 cm
525 cc	350-1485	450 cc	525 cc	15.0 cm	4.4 cm	630 cc	14.4 cm	5.6 cm
575 cc	350-1490	490 cc	575 cc	14.7 cm	4.8 cm	690 cc	14.2 cm	6.0 cm

Siltex® Round Spectrum®

MINIMUM FILL

MAXIMUM FILL

Device Volume	Catalog Number	Temp. Min. Volume	Final Min. Volume	Min. Diameter	Min. Projection	Final Max. Volume	Max. Diameter	Max. Projection
125 cc	354-2410M	105 cc	125 cc	9.4 cm	3.0 cm	150 cc	9.2 cm	3.5 cm
175 cc	354-2420M	150 cc	175 cc	10.5 cm	3.1 cm	210 cc	10.4 cm	3.9 cm
225 cc	354-2430M	190 cc	225 cc	11.3 cm	3.3 cm	270 cc	11.1 cm	4.2 cm
275 cc	354-2440M	235 cc	275 cc	12.0 cm	3.6 cm	330 cc	11.4 cm	4.4 cm
325 cc	354-2450M	275 cc	325 cc	12.9 cm	3.7 cm	390 cc	12.6 cm	4.5 cm
375 cc	354-2460M	320 cc	375 cc	13.4 cm	3.9 cm	450 cc	13.1 cm	4.8 cm
425 cc	354-2470M	360 cc	425 cc	14.0 cm	4.2 cm	510 cc	13.6 cm	5.0 cm
475 cc	354-2480M	400 cc	475 cc	14.5 cm	4.5 cm	570 cc	14.2 cm	5.7 cm

Note: Individual implant dimensions may vary slightly in products of this type. Not all units will conform exactly to the dimensions noted above.

Siltex® Contour Profile® Spectrum®

MINIMUM FILL

MAXIMUM FILL

Nominal Volume	Catalog Number	Temp. Min. Volume	Final Min. Volume	Min. Width	Min. Height	Min. Projection	Final Max. Volume	Max. Width	Max. Height	Max. Projection
275 cc	354-2511	235 cc	275 cc	11.5 cm	9.5 cm	5.1 cm	330 cc	11.1 cm	9.4 cm	6.3 cm
350 cc	354-2512	300 cc	350 cc	12.3 cm	10.5 cm	5.3 cm	420 cc	11.9 cm	10.1 cm	6.9 cm
450 cc	354-2513	380 cc	450 cc	13.2 cm	11.0 cm	6.1 cm	540 cc	12.7 cm	11.0 cm	7.5 cm
550 cc	354-2514	470 cc	550 cc	14.0 cm	11.9 cm	6.4 cm	660 cc	13.7 cm	11.9 cm	7.9 cm
650 cc	354-2515	550 cc	650 cc	15.0 cm	12.7 cm	6.6 cm	780 cc	14.5 cm	12.5 cm	8.2 cm

INJECTION DOMES

Micro Dome 350-MICRO

Includes stainless steel connector

Standard Dome 350-STD

Includes stainless steel connector

Dome Pack 350-DOMPK


Includes Standard and Micro Domes and both stainless steel and True Lock® connectors



Standard Dome





Micro Dome



Stainless Steel Connector



True-Lock® Connector System



201 Mentor Drive
Santa Barbara, CA 93111 USA
(805) 879-6000

Customer Service:
(800) 235-5731
Fax: (805) 967-7108

www.mentorcorp.com

Manufacturer:
Mentor
3041 Skyway Circle North
Irving, TX 75038 USA

European Representative:
Mentor Medical Systems B.V.
2333 CL Leiden, The Netherlands

INTERNATIONAL OPERATIONS

Mentor Medical Systems, Australia, Pty. Ltd. Tel: +61 2 9251 8111 Toll Free: +1 800 338 160 Fax: +61 2 9251 8333

Mentor Medical Systems, Canada, Inc. Tel: +1 905 725 7763 Toll Free: +1 800 668 6069 Fax: +1 905 725 7340

Mentor Deutschland GmbH Tel: +49 811 600 50 0 Fax: +49 811 600 50 11

Mentor Medical Systems France, S.A. Tel: +33 1 30 07 07 60 Fax: +33 1 30 07 07 67

Mentor Medical Systems Iberica, S.L. Tel: +34 91 562 2700 Fax: +34 91 728 3367

Mentor Medical Italia, S.r.l. Tel: +39 0 2 880 7761 Fax: +39 0 2 7209 3710

Mentor Medical Systems Ltd., U.K. Tel: +44 1635 511800 Fax: +44 1635 511810